

Ovum pick-up

A BIMONOCYTE B BIMONOCYTE-N

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Single lumen needle for oocyte aspiration and follicle flushing

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For atraumatic oocyte aspiration in IVF procedures. Single lumen, double-way needle, one way for oocyte aspiration (5), second way for follicle flushing (3).

Echogenic tip with atraumatic bevel (1). Oocyte aspiration is performed through a blue-colored tubing set that protects the oocyte from UV rays. The aspiration tubing set ends into a thermoplastic rubber bung (6) that is securely connected to the test tube. The needle is connected to both the aspiration and flushing tubings through a Y connector (2) (serving as finger grip) provided with two universal adapters that make the needle extremely versatile. Follicle flushing is performed through a plastic tubing set (3) ending in a stopcock (4) preventing any possible fluid reflux. Oocyte aspiration is controlled by a tubing clamp (10). Two different connections to the suction pump are available:

- short connection in **BIMONOCYTE (A)**, through a luer lock adapter (7);
 - long connection in **BIMONOCYTE-N (B)**, through an extension tubing 80 cm long, ending in a luer-lock adapter (8); both versions are supplied with an additional universal adapter (9) allowing multiple connections. Beta irradiated to eliminate any risk of toxicity from residual contamination common with EtO.
- BIMONOCYTE and BIMONOCYTE-N are MEA (mouse embryo assay) and LAL (endotoxin) tested.**

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10 pcs. per box. Single packaged, single use, Beta sterilized and CE marked. It is available with a 100 cm extension tubing upon request.

| Cod. A | Cod. B | Misura e lunghezza ago Needle size and length |
|--------|-----------|--|
| 770304 | 770304/TL | 16 G x 330 mm |
| 770302 | 770302/TL | 17 G x 330 mm |